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1. HB&L™ CULTURE KIT
2. HB&L™ P.A.R. TEST KIT
3. HB&L™ DEB KIT (Difficult Element Broth)

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


**Sterile**

CE Marked


CSF  
 Synovial Fluid  
 Pleural Fluid  
 Ascitic Fluid  
 Peritoneal Fluid


Synovial Fluid

CSF

Peritoneal Fluid


Pleural Fluid

Ascitic Fluid

Increased amount of fluid between abdominal structures

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
1. POSITIVE SAMPLE SELECTION AND MICROSCOPE EVALUATION

2. DILUTION 1:10 OF POSITIVE BLOOD CULTURE BROTH IN SALINE SOLUTION

3. 500 µl SUSPENSION DILUTION IN A HB&L™ CULTURE VIAL + DEB SUPPLEMENT UNTIL THE 0.5 MCFARLAND

4. ALIFAX AST with a customized antibiotic panel following:

- Gram indication
- Administered drug information
- Therapeutic protocols data and guide line

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From 1992 more than 60 scientific publications on Alifax technology.

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PERFORMANCE OF URO-QUICK, A NEW AUTOMATED  
METHOD FOR THE DETECTION OF BACTERIURIA

O. Soro, M. Raggi, G.C. Schito, Institute of Microbiology, University of Genoa - Italy

**EVALUATION OF AUTOMATED BACTERIURIA SCREENING SYSTEM  
IN SAMPLES COLLECTED IN THE PRESENCE OF BACTERIOSTATIC SUBSTANCES**

*I. Russo, E.M. Magliano, Microbiology Laboratory,  
NIGUARDA HOSPITAL, Milan, Italy*

## Urine culture automation new diagnostic and organizational paths

Presentation of Dr. Ricci L.

A.O. S.M. Nuova Microbiology Laboratory of Reggio Emilia (Italy)

# ANTIBIOGRAMMA

## Antibiotic Susceptibility Tests Directly on Urine Samples Using Uro-Quick, A Rapid Automated System

S. ROVETA - A. MARCHESE - E.A. DEBBIA

Institute of Microbiology, C.A. Romanzi, DISCAT- University of Genoa - Largo Rosanna Benzi 10 -16132 Genoa, Italy.

- » **12.579 urine samples** collected from in-patients (9906) and out-patients (2673), June 2003 - March 2004. Positive samples were processed for antimicrobial susceptibility tests by **standard disk-diffusion method**.
- » Globally **1590 gram-negative** and **358 gram-positive** pathogens were isolated. Polymicrobial samples (3.5%) were excluded from the study.
- » **The present findings demonstrate and confirm, on a large number of bacteria, that Uro-Quick system performs efficiently in the determination of the antibiotic susceptibility patterns of Gram- and Gram+ bacteria representing the most important community acquired and nosocomial urinary tract pathogens with an agreement > 90% with the standard disk-diffusion method”.**

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Published: 2009.02.01

## Authors' Contribution:

- A** Study Design
- B** Data Collection
- C** Statistical Analysis
- D** Data Interpretation
- E** Manuscript Preparation
- F** Literature Search
- G** Funds Collection

## A novel culturing system for fluid samples

Carla Fontana<sup>1,2A,B,C,D,E,F,G</sup>, Marco Favaro<sup>1B,C,D,F</sup>, Silvia Minelli<sup>2B,D</sup>, Maria C. Bossa<sup>2B,D</sup>,  
Anna Altieri<sup>2B,D</sup>, Cartesio Favalli<sup>1,2G</sup>

<sup>1</sup> Department of Experimental Medicine and Biochemical Sciences, "Tor Vergata", University of Rome, Rome, Italy

<sup>2</sup> Clinical Microbiology Laboratories, Polyclinic of Tor Vergata, Rome, Italy

**Source of support:** Departmental sources

# ANTIBIOGRAMMA IN EMOCOLTURE POSITIVE

MICROBIOLOGIA MEDICA, Vol. 25 (1), 2010

RAPID ANTIBIOTIC SUSCEPTIBILITY TESTS

## Preliminary indications for antibiotic susceptibility tests in less than six hours in positive blood cultures

**Vesselina Kroumova, Elisa Gobbato, Paola Macaluso, Sabrina Tamburelli, Francesca Marini, Marcella Perone, Stefania Orlandi, Morena Viviani, Giacomo Fortina**

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**European Society of Intensive Care Medicine (ESICM)  
Annual Congress**

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MICROBIOLOGIA MEDICA, Vol. 25 (2), 2010

OSTEOARTICULAR SAMPLE CULTURE IN LIQUID BROTH (HB&L SYSTEM)

## Evaluation of the HB&L system for the culture of prosthetic and osteoarticular origin samples

**Agostina Ronca<sup>1</sup>, Sabrina Brenci<sup>1</sup>, Giuliana Carrega<sup>2</sup>, Gianni Riccio<sup>2</sup>, Luisa Santoriello<sup>1</sup>**

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**Abstract: O538**

Citation: Clinical Microbiology and Infection. Volume 16 Supplement No. 2, Page S120

## **Speeding up identification of bacteria from urine by the use of an Alifax HB&L Uroquattro incubator followed by Maldi Biotyper identification**

M. Kostrzewa, C. Boogen, J. Maier, U. Weller

*(Bremen, Cologne, DE)*

**Conclusion:** This protocol enables reliable bacterial identification together with a bacterial count in few more than 3 hours in cases of mono-bacterial infection using the standard protocol of the HB&L. Ongoing improvement to the MBT identification algorithms will resolve the problem of poly-bacterial infections. Further extension of this approach to usually sterile body fluids, tracheal secretions and pleural aspirates also here may lead to a dramatic reduction in turnaround time of microbiological labs for identification and subsequent AST-testing.

**ORAL COMMUNICATION ECCMID 2010** 

## Rapid Identification of Bacteria from Urine using the ALIFAX HB&L Uroquattro incubator followed by MALDI Biotyper identification

<sup>2</sup>Markus Kostrzewa, <sup>1</sup>Christiane Boogen, <sup>1</sup>Julia Maier, and <sup>1</sup>Ulrich Weller

<sup>1</sup>Labor Dr. Boogen, Hauptstr. 71-73; Cologne, Germany  
<sup>2</sup>Bruker Daltonik GmbH, Fahrenheitstr. 4, Bremen, Germany

### Introduction

Direct identification of bacteria from infected urine by the MALDI Biotyper has been shown previously. To enable quantification as well as identification, which is necessary for routine diagnostic, we have tested the combination of an automated quantification system, the Uroquattro (ALIFAX HB&L), with the MALDI-TOF workflow.

Fig.1: General workflow of MALDI-TOF MS microbial identification

### Methods

322 urine samples from routine diagnostic of out patients were analyzed in parallel.

*Standard workflow in Laboratory Dr. Boogen:*

- Quantitative culture, over night (16h)
- Potentially isolation of pure colonies for 6-16h
- Identification with MALDI Biotyper (minutes)
- Merlin Micronaut AST (6-16h)

*Workflow for the proof of concept study:*

A 500 µl aliquot was transferred into a 2ml HB&L tube with nutrient broth and stirrer, placed into the Uroquattro instrument for 3h at 37°C and the CFU value was read. Thereby, the quantification of bacteria was performed and, in parallel, the bacteria were amplified to a number which is well suitable for MALDI-TOF analysis. After incubation, 1ml aliquots from the positive vials were harvested by centrifugation (5 min at 15000 g in a 1.5 ml Eppendorf cup). Supernatant was discarded and the pellet resuspended in 1 ml deionized

water and re-centrifuged in the same manner. The supernatant was discarded again and part of pellet transferred onto a MALDI target, air dried, overlaid with matrix solution (HCCA), again air dried and subjected to identification by mass spectrometry. Profile mass spectra were acquired in the linear positive mode (2000-20000 Dalton) and subjected to analysis with the MALDI Biotyper 2.0 software (database with 1857 species).


Fig. 2: Novel workflow combining the Uroquattro system with the MALDI Biotyper identification.

### Results

In a prospective proof-of-concept study with 322 urines 145 samples were found to be culture-positive over a period of 4 weeks. These samples were found to have a bacterial count over 10<sup>4</sup> by both the standard workflow of the laboratory and the automated quantification in the Uroquattro system.

Monoinfections were found in 127 cases (87.5%), both by the standard method and the new workflow with total agreement of results. The following species were detected:

*E. coli* (67), *P. mirabilis* (12), *E. faecalis* (8), *K. pneumoniae* (8), *P. aeruginosa* (6), *R. ornithinolytica* (6), *M. morgani* (6), *S. aureus* (4), *E. aerogenes* (3), *C. koseri* (2), *K. oxytoca* (2), *S. marcescens* (2), and *S. enterica* (1)

18 samples (12,5%) were found to be infected by more than one microorganism. In these cases the novel method did not come to a conclusive result (no reliable identification score).

ASM 2010, San Diego  
C-1587

### Conclusions:

More than 80% of samples in Laboratory Boogen are mono-infections. Currently, more than 90% of these infections can be identified by the HB&L Uroquattro and MALDI Biotyper.

Following growth in the HB&L and identification by the MALDI Biotyper, the residual samples are used for AST (Merlin Micronaut), McFarland continuously measured by the HB&L.

The whole procedure can be completed in a working day.

### Perspectives

MALDI Biotyper and Uroquattro:

- Automated analysis of mixed cultures is being developed
- Extension of procedure to primary sterile body fluids, tracheal secretions and pleural aspirates

MALDI Biotyper

## Assessment of a light-scattering system in the culture screening of biological fluid samples

**Paolo Lanzafame**<sup>1</sup>, **Maira Zoppelletto**<sup>2</sup>, **Marina Gaino**<sup>1</sup>, **Patrizia Ober**<sup>1</sup>

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2011 Trends in microbiology

### Summary

The diagnosis of invasive infections is commonly based on the cultures of pathogenic microorganisms; the morphological and biochemical identification is followed by several types of antimicrobial susceptibility tests.

Early identification and rapid antimicrobial susceptibility testing of microorganisms, particularly in the case of serious infections, represent high priorities in clinical microbiology. Our study aimed at evaluating a semiautomated system (HB&L by Alifax S.p.A., Padova, Italy) for the recovery of microorganisms from sterile body fluids other than traditional cultural methods.

HB&L offers the possibility for the culture of an additional 5.8% samples, if compared with the reference cultures. The two methods showed good concordance (85%), with high sensitivity (97.87%) and specificity (86.3%) and very good negative predictive value (94.83%) and a positive predictive value (82.14%).

The microbial counts showed excellent agreement in the results between HB&L and traditional methods.

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Tabella 2 Alfred60 vs Piastra

Table with 2 columns: Metric and Value. Rows include SENSIBILITA' (99.8%), SPECIFICITA' (90.0%), VPN (99.9%), VPP (83.6%), and Concordanza (93.3%).

Tabella 3 Alfred60 vs Gold Standard

Table with 2 columns: Metric and Value. Rows include SENSIBILITA' (99.8%), SPECIFICITA' (87.1%), VPN (99.9%), VPP (78.1%), and Concordanza (91.1%).

Tabella 4 Piastra vs Gold Standard

Table with 2 columns: Metric and Value. Rows include SENSIBILITA' (96.2%), SPECIFICITA' (95.0%), VPN (98.2%), VPP (89.9%), and Concordanza (94.1%).

Received: 29 March 2011

Revised: 26 May 2011

Accepted: 26 May 2011

Published online in Wiley Online Library

*Rapid Commun. Mass Spectrom.* 2011, 25, 2247–2249  
(wileyonlinelibrary.com) DOI: 10.1002/rcm.5113

## Direct identification of bacteria in blood culture by matrix-assisted laser desorption/ionization time-of-flight mass spectrometry: a new methodological approach

Vesselina Kroumova<sup>1\*</sup>, Elisa Gobbato<sup>1</sup>, Elisa Basso<sup>2</sup>, Luca Mucedola<sup>1</sup>, Tommaso Giani<sup>3</sup>  
and Giacomo Fortina<sup>1</sup>

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<sup>2</sup>Bruker Daltonik, S.r.l., Via Cluentina 26/r, Macerata 62100, Italy

<sup>3</sup>Department of Biotechnology, Section of Microbiology, University of Siena, Viale Bracci 2, Siena 53100, Italy

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

2012 SIM e ISF

### From positive blood culture to microbiological diagnosis in four hours by MALDI-TOF mass spectrometry bacterial identification and rapid antibiogram

S. Barnini, V. Brucculeri, P. Morici, A. Lupetti, M. Campa  
Azienda Ospedaliero-Universitaria Pisana, Microbiology Unit, Pisa, Italy

### Background

Time is precious for patients and clinicians facing septic events. In the recent past years, many efforts have been done to develop rapid and reliable tests, especially by molecular methods, which still have several limits in determining antimicrobial susceptibility. Microbiological diagnosis of sepsis is still based on blood culture, which allows isolation, identification and determination of antibiotic susceptibility. Nowadays, MALDI-TOF mass spectrometry (MS) allows almost instantaneous bacterial identification, even directly, from positive blood cultures, but antibiotic susceptibility results are usually delayed by 24-48 hours, compared to the timing of microorganism identification. Since inadequate initial antimicrobial treatment is associated with worse outcomes for patients, in this study attention was paid on tuning a rapid susceptibility test, which could give reliable answers in few hours, after bacterial identification obtained from positive blood cultures by means of two short methods applied to MALDI-TOF MS.

# PUNTE DI CATETERE VENOSO CENTRALE

## Improved diagnosis of central venous catheter-related bloodstream infections using the HB&L UROQUATTRO™ system

C. Fontana • M. Favaro • M. C. Bossa • S. Minelli •  
A. Altieri • M. Pelliccioni • F. Falcione • L. Di Traglia •  
O. Cicchetti • C. Favalli

[Eur J Clin Microbiol Infect Dis.](#) 2012 Nov;31(11):3139-44.  
doi: 10.1007/s10096-012-1676-9. Epub 2012 Jun 27

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